







**Secure Service**

Rock-solid security is provided for valuable patient data, without requiring firewall modifications and without requiring a VPN.

APPROACH	VALUE
Secure connectivity via a web service	Innovative solution eliminates the need for cumbersome VPNs
128- to 256-bit SSL encryption	Layered security ensures protection of confidential patient data
Connectivity to remote sites and EMR systems without requiring modifications to the firewall	No adjustments to existing security configurations are needed

**Scalable Solution**

Built on an industry-proven integration platform, Corepoint Community Exchange expands your relationships in the healthcare community without sacrificing performance or increasing operational costs.

APPROACH	VALUE
Leverages your Corepoint Integration Engine investment	Reduces the cost of deployment and enables robust interfacing through a trusted integration platform
Uses a "locking" mechanism as a security check when connectivity is scaled	Maintains First In First Out (FIFO) order in each message stream
Flexible to work with all data types	Works with both HL7 and XML messages, and populates a CCR or CCD as required
Builds on existing connectivity	Utilizes existing configurations to quickly expand the number of connections with minimal modifications at the central server
High volume	Can connect to hundreds of remote sites or EMR systems without sacrificing performance